

TDTR_vH3: rapid data processing,
analysis, and note-keeping for
general thermal modeling.

by Greg Hohensee, based on Joe
Feser's MATLAB scripts and beam-
offset code.


Why?

- In weekly meetings, I was expected to:
 - Re-analyze whole data sets (sometimes!)
 - Not mix up thermal parameters by accident.
 - Fit data manually for intuition.
 - Show ALL thermal model parameters and plot fit results in a specific format.
 - For low temperature measurements, self-consistently update thermal parameters for steady-state heating.

Features Overview

- Two front-end scripts for each TDTR session: your analysis saved, no back-end customization required.
- Assisted phase & time shifts, ps acoustics and Brillouin reading.
- Integrated database for your materials' thermal parameters.
- Manual fitting to any signal, with option for sensitivity calculations.
- Self-consistent steady-state temperature adjustment.
- Generates and saves detailed fit result figure.
- NEW: Parametric sensitivity plots

Script Hierarchy



Parameter Bundles, or: I should have made objects instead of cell arrays.


- datparams - {tdelay data datadir offset}
- matparams - {LCTE aniso BI n_toplayer TCR doughnut}
- sysparams - {tau_rep f r_pump r_probe}
- calparams - {Zind sigfit intscheme nnodes consider_error LCTE_err T0_err P0_err}
- Tparams - {T0, T_LCTE, A_pump, A_probe, absC, perpulse, jabs, jtrans}
- INITIALIZE_CELLPARAMS_vH3.m contains extensive comments explaining all these parameters. It also handles unpacking and repacking of these cell parameter arrays for many other scripts; update this script if you make changes to the cell arrays.

Front-end: Process and Analyze

- Easy documentation and revisiting: just save new (dated) copies of `process_yymmdd.m` and `analyze_yymmdd.m` scripts for each TDTR session.
- **Process** helps setting zero of delay time, setting the phase, picking picosecond acoustics, and examining the signals of your data series at particular delay times.
 - Saves phase- and time-shifted data with the “_shifted.txt” suffix to “yymmdd_edit” folder.
- **Analyze** lets you set all system and thermal parameters, and customize many aspects of the modeling and fitting process. Model parameters are specified as arrays, so it can loop through a full TDTR session and save the solutions.

SetPhase.m

- Shows short- and long-time delay $V_{\text{out}}(t)$, raw versus shifted.
- Options:
 - Manual (Joe Feser's code)
 - Automatic (compare average $V_{\text{out}}(t)$ before/after 0 ps)
 - Linear fit (for seeing ratio without $V_{\text{out}}(t)$ noise)
- Other process scripts:
 - **set_t0.m** – helps you set zero of time delay.
 - **PickEchoes.m** – lets you pick picosecond acoustic timing.
 - **PickEchoesB.m** – helps with Brillouin signals.
 - **AutoSetPhase.m** – Auto shifts data without user input.



TDTR_MANUALFIT_vH3.m: manual fitting for intuition

```
Pressure is 3.9 GPa  
Manual fitting to ratio -V(in)/V(out)...
```

```
Z =  
  
    0.4373
```

```
Mat =
```

```
Current material fit parameters LCTE(i,j):  
L(2) = 0.1000 W/m-K  
L(1) = 7.0000 W/m-K
```

```
Enter 1 if done, 2 for a sensitivity plot, 3 to rescale; else hit "Enter": 4  
Hey! Invalid input. Go home, you are drunk.  
Adjust parameter L(2): 0.15  
Adjust parameter L(1): 9
```


<input type="checkbox"/>	PP
<input type="checkbox"/>	T
<input type="checkbox"/>	T0
<input type="checkbox"/>	TCR
<input type="checkbox"/>	TDebye
<input checked="" type="checkbox"/>	T_LCTE
<input checked="" type="checkbox"/>	Tparams
<input type="checkbox"/>	VinVout
<input type="checkbox"/>	Vin_data
<input type="checkbox"/>	Vin_raw
<input type="checkbox"/>	Vout_data
<input type="checkbox"/>	Vout_raw
<input type="checkbox"/>	Voutlinfit
<input type="checkbox"/>	XPsubs
<input type="checkbox"/>	Xguess

- Can generate a sensitivity plot at any time.
 - To speed up sensitivity plotting, go into senseplot_offset_vH3.m script and edit “LCTE_sens_consider” and “sys_consider” variables to select which sensitivities to calculate, and which ones to skip.
- “3 to rescale”: manually shift normalization reference level for $V_{in}(t)$ and $V_{out}(t)$ fitting, in case of noise or acoustics.

Sensitivity conventions and coupled thermal parameters

- Sensitivity (*and error bars*) are **never** computed relative to raw data.
 - Sensitivities and error bar computations use perturbations from your best-fit model: this avoids extrinsic effects from noise and second-guessing your manual fitting judgment.
- Anisotropy “eta” = Λ_x / Λ_z , in-plane divided by cross-plane thermal conductivity. Λ_z is the “lambda” variable that you’re familiar with.
 - Convention: if model layer is anisotropic, perturbing Λ_z also updates eta so that Λ_x remains fixed relative to Λ_z . Sensitivity to eta, $S(\text{eta})$, is equivalent to $S(\Lambda_x)$.
- Transducer absorption layers:
 - Program detects absorption and transducer layers (user specifies this in analyze script!). Perturbations are done for the transducer as a whole, not the two model layers separately.

Temperature-dependent thermal parameters, self-consistently




TDTR_TTM_vH2: Two-channel model for *substrates only*

- Two-channel thermal model designed by Rich.
 - See: R. B. Wilson *et al.*, *Phys. Rev. B.* **88**, 144305 (2013).
- Can fit thermal model to the ratio of two TDTR measurements at two different modulation frequencies. Also sensitivities.
 - See: G. T. Hohensee *et al.*, *Phys. Rev. B.* **89**, 024422 (2014).
- Designed for transducer on bulk spin-ladder sample: for multilayer 2- or 3-channel thermal modeling, ask Rich.
- User may need to customize anisotropy in two-channel substrate using the back-end scripts.
- TTM_vH2 has *not* been updated to vH3.

parametric_senseplot_vH3.m:

Parametric Sensitivity Plotter

- Controlled by `sense_setup.m`, or `sense_setup_retro.m` for users of Joe Feser's TDTR_MAIN script.
- $S(X) = d\log(R)/d\log(X)$ versus model parameter Y or offset, not just time delay.
- R can be ratio, normalized $V(\text{in})$ or $V(\text{out})$, beam offset, or FWHM of beam offset.
- Verified against Cahill group publications
 - David Cahill *et al.*, "Nanoscale Thermal Transport II," *Appl. Phys. Rev.* **1**, 011305 (2014)
 - Yee Kan Koh, *Ph.D. Thesis*, UIUC (2010)
 - Xiaojia Wang *et al.*, *Phys. Rev. B.* **88**, 075310 (2013)
 - But not beam offset! Proceed with caution, verify against published sensitivities [Such as: J. Feser *et al.*, *Rev. Sci. Instrum.* **85**, 104903 (2014)]



Thank you for your attention!

Slides are be posted on GitHub along with the TDTR_vH3 package.

http://github.com/gthohensee/TDTR_vH3_pub